

PATENT APPLICATION

INVENTORS:	Name	Residence
	James D. Browning	Columbus, Ohio
	Vincent D. McGinniss	Sunbury, Ohio
	Steven M. Risser	Reynoldsburg, Ohio
	Dale L. Robinson	Powell, Ohio
	Robert E. Russell, III	Newark, Ohio
	Robert Boyd	Dublin, OH
	Matthew Fleming	Upper Arlington, OH

TITLE: SKIN CONDUCTION AND TRANSPORT SYSTEMS

ATTORNEY Klaus H. Wiesmann

REG. NO.: 30,437

CORRESPONDENCE
ADDRESS: Battelle Memorial Institute
505 King Avenue
Columbus, Ohio 43201-2693
(614) 424-7798 or 6589

The invention was not made by an agency of the United States Government nor under contract with an agency of the United States Government.

Express Mail No. EU417945385US

SKIN CONDUCTION AND TRANSPORT SYSTEMS

5

This application claims benefit of provisional application 60/393,036 filed June 28, 2002.

The entire disclosure of the referenced provisional application is incorporated herein by reference as if completely rewritten herein.

10

FIELD OF THE INVENTION

This invention relates to materials and methods for providing enhanced mucous membrane or skin conductivity in electrical measurements involving the surface of an animal or human. The invention also provides materials and methods for enhanced transport of substances across the mucous membrane or skin boundary of an animal or human.

BACKGROUND OF THE INVENTION

The measurement of internal body conditions by means of external measurements is the preferred measurement technique for diagnosing disease or monitoring body conditions. The goal for measurement by external measurements is to be as minimally invasive as possible to achieve measurement of the desired internal condition. Many devices have been designed to achieve the above goals through the use of electrical or electrochemical surface measurements.

Another important technique that is becoming increasingly important is the transdermal administration of drugs as a preferred mode of delivery. Transdermal delivery of drugs provides many advantages over other means of administration such as oral or by injection. Advantages include safety, convenience, increased patient compliance, non-interrupted therapy, avoidance of the hepatic first pass effect, and the high degree of control of any administered drug.

Many drugs are not amenable to transdermal delivery due to the well known barrier properties of the skin. Drug molecules that are intended to penetrate into the body through intact skin must first penetrate the stratum corneum and materials therein or on its surface. The drug molecule must
5 then penetrate viable epidermis, the papillary dermis, and then the capillary walls before entering systemic circulation. In its path into the body, each of the mentioned tissues will exhibit different resistances to penetration by the same drug molecule. It is the stratum corneum, however, that presents the greatest barrier to absorption by topical administration or transdermal
10 administration for most molecules. Mucous membranes present similar problems depending on the materials that are desired to penetrate the membrane.

BRIEF DESCRIPTION OF THE INVENTION

15 Broadly, one embodiment of the invention provides for enhanced electrical contact between an electrode and a patients surface (e.g. mucous membrane, skin). Enhanced electrical contact typically provides for greater reliability of measurements for the diagnosis of body conditions. In addition, improved electrical contact typically provides for greater sensitivity, comfort
20 and ease of use for many applications.

Broadly, another embodiment of the invention provides for enhanced transdermal or trans-mucous membrane delivery of drugs or other materials.

Broadly, a still further embodiment of the invention provides for diagnosis of bodily conditions by transdermal or trans-mucous membrane
25 movement of body fluids to an external detection device.

A first broad embodiment for an electrode providing electrical contact with a surface of a patient includes a conductive member; and a conduction enhancer in contact with the conductive member including a carrier and a protein/fatty acid based compound. Typically, the protein/fatty acid based
30 compound may be a lipopolypeptide, the lipopolypeptide may be an acyl peptide, and the acyl peptide may be a material selected from the group consisting of Lamepon S™, MayTein C™, MayTein CT™, and mixtures thereof.

The carrier and conduction enhancer typically provide electrical contact with an electrical resistivity of less than about 10K Ohms when the electrode is applied to the patient's surface.

Typical activities for the surfactant (conduction enhancer) according to the invention include: activity between about 0.25% and about 60%; preferably an activity between about 4% and about 50%, more preferably an activity between about 5% and about 30%; and most preferably has an activity between about 10% and about 30%.

A yet further embodiment of the invention includes a method for decreasing the electrical resistivity between an electrode and the surface of a patient including placing a carrier and a conduction enhancer of a protein/fatty acid based compound between the electrode and the surface. Typically, the method according the carrier and conduction enhancer provide electrical contact with an electrical resistivity comprising less than about 10K Ohms when the electrode is applied to the patient's surface. Typically, the electrical resistivity is obtained in about 0.001 seconds to about 3 minutes; in other embodiments the electrical resistivity is obtained within about 0.01 seconds to about 30 seconds. In yet other embodiments the electrical resistivity is less than about 6K Ohms. Typically, the electrical resistivity is maintained for at least about 8 hours, and most preferably the electrical resistivity is maintained for at least about 72 hours. Typically, the composition comprises a gelling agent.

In some embodiments, the invention provides for a composition for enhancing the electrical conductivity between an electrode and a patient's surface including: a mixture comprising a carrier and a protein/fatty acid based compound. The composition may be a lipopolypeptide, an acyl peptide and the like. Typically, the acyl peptide comprises a material selected from the group consisting of Lamepon S™, MayTein C™, MayTein CT™, and mixtures thereof. Typically, the carrier and conduction enhancer provide electrical contact with an electrical resistivity comprising less than about 10K Ohms when the electrode is applied to the patient's surface.

In further embodiments, the invention provides for an electrode providing electrical contact with a surface of a patient, the electrode comprising: a conductive member; and a conduction enhancer in contact with the conductive member comprising a carrier and a surfactant represented by
5 the formula:



wherein R, R', and R'' are the same or different and may be independently
10 selected from the group consisting of alkyl, aryl, amine, carbonyl, and carboxyl moieties; R, and R'' may also be independently selected from the group consisting of -H, and -SH;
wherein the repeat unit, n, is an integer from about 2 to about 2000; and wherein M is a metal ion. Typically, when a carbon containing moiety is
15 selected, R, R', and R'' have 1-20 carbon atoms. Typically, the repeat unit, n, is an integer from about 150 to about 1800. In some embodiments, the surfactant is a mixture of compounds selected from the formula. The metal ion, M, is typically selected from the group consisting of K⁺, Na⁺, and mixtures thereof.

20 In additional embodiments, the invention includes a method for making electrical contact between an electrode and a patient's surface, which comprises the step of applying an electrode having a surface coated with a mixture of a carrier and a conduction enhancer comprising a protein/fatty acid based compound. Typically, the patient's skin is unabraded. Some
25 embodiments provide for conditions where the electrical resistivity is obtained in about 0.001 seconds to about 3 minutes, more preferably where the electrical resistivity is obtained within about 0.01 seconds to about 30 seconds. Most preferably the electrical resistivity is less than about 6K Ohms. The resistivity is preferably maintained for at least about 8 hours and most
30 preferably at least about 72 hours.

In another embodiment according to the invention, there is provided a method for the noninvasive measurement of body substances from a patient

including the steps of applying a multilayer patch to the skin of the patient; wherein at least one layer of the multilayer patch comprises a transdermal migration-enhancing amount of an acyl peptide; and measuring a body substance, the transdermal migration of which into the multilayer patch is facilitated by the acyl peptide. The method may include the additional step of enhancing the migration of the body substance by reverse electro transport. In some embodiments, the reverse electro-transport comprises iontophoresis.

In yet another embodiment of the invention, there is provided a device for measuring the quantity of a body substance in a patient comprising: a multilayer patch device including a first layer comprising an acyl peptide; a second layer comprising a reaction layer; and a third layer comprising a readout visible to a user indicative of the level of the body substance in the patient.

A further embodiment includes a method for enhancing patient surface/electrode conduction for electrosurgery in a patient comprising: applying a conduction enhancing amount of a mixture of a carrier material and a conduction enhancer comprising a protein/fatty acid based compound.

An additional embodiment provides for a surgical electrode for electrosurgery comprising a conductor and a conduction-enhancing amount of an acyl peptide. And another embodiment provides for an electrode for electrosurgery in a patient comprising: a conductor; and a conduction-enhancing amount of a mixture of a carrier and a surfactant selected from the group consisting of Lamepon S™, MayTein C™, MayTein CT™, and mixtures thereof, wherein the mixture is applied to a surface of the conductor.

25

BRIEF DESCRIPTION OF THE DRAWINGS

Figure 1 is a schematic drawing of a top view of apparatus used for testing surfactant penetration of surfaces such as membranes or skin.

Figure 2 is a schematic drawing of a side cutaway view of the apparatus of Figure 1 along section 2.

Figure 3 is a schematic drawing of a side cutaway view of another embodiment of the apparatus for determining conductivity.

Figure 4 is a schematic drawing of a multilayer patch according to the invention.

Figure 5 is a schematic drawing of an electrosurgery electrode according to the invention.

5 DETAILED DESCRIPTION OF THE INVENTION AND BEST MODE

Definitions

% Active - this term means that the stated percent (weight percent) of the substance is contained in the solution or gel.

10 % - all percent measurements are in weight percent unless stated otherwise.

Breakthrough time - the time required for the solution to be transported through the membrane(s), and contact and bridge two pieces of copper wire (wire gap ~ 0.0625 inch).

15

Conductive – as used herein refers to electrical conductivity unless otherwise explicitly stated.

Non-conductive – as used herein refers to electrically non-conductivity unless otherwise explicitly stated.

20

Patient's surface – in one embodiment the patient's surface to which an electrode or a patch is applied is a mucous membrane while in another embodiment the patient's surface to which an electrode or patch is applied is the skin.

25

* * * * *

Broadly the invention discloses new materials and methods for the testing and evaluation of mammalian and particularly human body parameters. The parameters include electrical and chemical characteristics of the body. The present invention provides for improved conductivity between electrodes and a patient's (animal or human) mucous membrane or skin.

30

The materials of the present invention can be applied as liquids or gels (mixed with gelling agents). The liquid form is typically applied by absorption into foams, fibrous networks or other materials that absorb, adsorb or otherwise hold the material in place.

5 Typical materials useful in the invention include surfactants to enhance conductivity at the patient's surface/electrode interface. In other embodiments the materials of the invention provide enhanced drug delivery through the intact skin or mucous membranes with and without the aid of electrical current. Typically, the surfactant is a protein based surfactant
10 having hydrophilic and hydrophobic properties, and/or lipophilic and lipophobic properties. A typical protein based surfactant is a lipopolypeptide. The lipopolypeptide is typically made from long chain fatty acids and hydrolyzed protein. The long chain fatty acids may be of animal or plant origin. The polypeptide may likewise be of animal or plant origin. Preferably
15 the lipopolypeptide is a natural or modified coconut oil and hydrolyzed protein (typically collagen). The typical general chemical structure can be represented as follows:



20 Typically, R, R', and R'' are the same or different and may be selected from the group consisting of alkyl, aryl, amine, carbonyl, or carboxyl moieties; and R and R' may also be independently selected from the group consisting of -H, or -SH. In some embodiments, when a carbon containing moiety is selected,
25 R, R', and R'' typically have about 1 to about 20 carbon atoms; preferably R'' typically has about 2 to about 10 carbon atoms. The value for the repeat unit, n, is typically about 2 to about 2000 and preferably about 150 to about 1800. In some embodiments, the lipopolysacharide represented by the formula is a complex mixture having various R, R', R'' groups, that typically
30 have various amounts of carbon atoms, and various values for the repeat unit, n. The ion -M is typically a metal ion (such as K⁺, or Na⁺, or mixtures thereof) that provides acceptable surfactant properties.

The surfactant materials are typically prepared by the reaction of the carboxyl group (-COOH) of a fatty acid with the amino group (-NH₂) of a polypeptide. The length of the fatty acid chain and the length of the polypeptide chain as well as the side groups are selected to provide characteristics for matching the physical and chemical properties of the surface of the patient and (e.g. skin or mucous membrane) and enhancing electrical conductivity between the patient's surface and an electrode.

In some embodiments, the surfactant materials are typically acylated polypeptide salts. The acylated polypeptide salts are typically represented by Lamepon STM, MayTein CTM and MayTein CTTM. MayTein CTM is similar to Lamepon STM. MayTein CTTM has an additional R'' group located at the asterisk in the above formula. M may be a K⁺, Na⁺, or similar ion. Lamepons and MayTeins may be obtained from Henkel Corp. They are typically condensation products of protein hydrolysates and vegetable fatty acids. Lamepon STM is a lipopolypeptide made from natural coconut oil and hydrolized collagen.

Typically, for improved electrical conductivity, the surfactants of the present invention are provided having activities of between about 0.25% to about 70%. Preferably typical activity ranges are between about 4% and 50%. More preferably are activities between about 5% and 30%. At the lower end of the activity scale the activities are most preferably above 10%. This is because the lowest electrical resistance for the tested substances has been observed in the neighborhood of about 20%.

Test Apparatus for Artificial Skin

Referring now to Figure 1, this figure illustrates a schematic diagram of a top view of testing apparatus 100. Apparatus 100 was used for testing surfactant penetration of several membrane types with several surfactants. A base 102 for testing apparatus 100 is typically constructed of poorly electrically conducting or electrically non-conducting material such as plastic and the like. The particular material used in the tests was Plexiglas and was about 0.125 inch thick. A groove 103 is cut into base 102 so that it went

completely across the base 102. Two pieces of 20 gauge copper wire 112, 114 are inserted into groove 103 and placed such that a gap 116 of about 1/16" is maintained between the copper wires 112, 114 at about the centerline 117 of the base 102.

5 A test material 108 is typically placed over the gap 116 so that it is centered thereon. A cover 104 of similar or like material as base 102 was placed over the test material. Cover 104 has a centrally placed opening 106.

 When used for tests, a test material was placed so it covered the centrally placed opening. One or more drops of surfactant material were
10 placed dropwise onto a test material its upper surface 109 and substantially in the area over gap 116.

 A Fluke ohmmeter 110 (e.g. Model 89 IV True RMS Multimeter) was connected to copper wires 112, 114 via leads 122, 124 for testing artificial skin such as one or more layers of membranes represented by test material
15 108.

 Typical test materials used with the test apparatus are layers of polycarbonate and nylon as more fully discussed in the examples below.

 Four identical 100g weighted vials 132 containing weights 134 were placed in the corners of the cover 104 to hold the cover 104 over the test
20 material 108. Figure 2 as a sectional view of Figure 1 along the plane represented by Section 2-2.

 Referring now to Figure 3, in some tests the four vials 132 for weighting were replaced by one large vial 132-1. Otherwise the materials were the same.

25

Example 1 (Conductivity and Resistivity of Surfactant Solutions)

 Tests were made of sixteen (16) surfactants. Fifteen of the surfactants had been diluted to an activity of 5% with deionized water. A deionized water sample EEG-16 was used as a control. One sample EEG-14, consisting
30 of 3ML-15638 oil/water emulsion was determined to be already very diluted and so was used, as received, at 80% activity. The surfactant solutions were all mixed for 30 minutes prior to electrical conductivity testing using a VWR

conductivity meter Cat# 23226-50, Probe Cat# 23226-524, range 0.1-200 u S/cm. Results are shown in Table 1.

Table 1. Conductivity and Resistivity of Standardized Surfactant Solutions

5

Sample	Surfactant	As Rec-eived (% Active)	Surfactant grams @ 5%	De-ionized Water (gram)	Surfactant solution total weight (gram)	% Active	Surfactant solution conductivity VWR (μS/cm)	Surfactant solution Resistivity (Ohm-meter)
EEG-1	Actrosol C-75™ (clear)	70%	1.07	13.93	15	5	8860	1.1
EEG-2	Actrosol C-85™ (disp.)	75%	1.00	14.00	15	5	5300	1.9
EEG-3	Actrosol OY-75™ (disp.)	75%	1.00	14.00	15	5	5230	1.9
EEG-4	Chembetaine C™ (disp.)	35%	2.14	12.86	15	5	5880	1.7
EEG5	Desonic DA-4™ (disp.)	100%	0.75	14.25	15	5	12.1	826.4
EEG-6	DeSulf SLES 301™ (clear)	25%	3.00	12.00	15	5	6880	1.5
EEG-7	Dodecyl Sulfate Sodium Salt (clear)	98%	0.77	14.23	15	5	6300	1.6
EEG-8	Fluorad FC-129™ (clear)	50%	1.50	13.50	15	5	4820	2.1
EEG-9	Lamepon S™ (clear)	32%	2.34	12.66	15	5	78.3	127.8
EEG-10	Lodyne S-103A™ (clear)	45%	1.67	13.33	15	5	4890	2
EEG-11	Lodyne S-222N™ (disp.)	85%	0.88	14.12	15	5	211	47.4
EEG-12	MayTein C™ (clear)	36%	2.08	12.92	15	5	9000	1.11
EEG-13	MayTein CT™ (clear)	38%	1.97	13.03	15	5	4600	2.2
EEG-14	3M L-15638™ (emulsion (o/w)	80%	15.00	0.00	15	80	58	172.4
EEG-15	T-MAZ 20 Peg™ Ester (clear)	100%	0.75	14.25	15	5	129.5	77.22
EEG-16	Deionized Water	--	0.00	15.00	15	--	9.43	1060.4
EEG-17	High Purity Water	--	0.00	15.00	15	--	1.51	6622.5

Example 2 - (Polycarbonate Skin Model)

In Example 2, tests were performed on surfactants shown in Table 2
 10 using test apparatus 100. The test material 108 consisted of five layer stack

of polycarbonate membranes. The membrane consisted of a polycarbonate Whatman Nuclepore™ having a 0.05 µm pore size (Cat. # 111703). Each had a thickness of 6 µm resulting in a stack about 30 µm thick. A Fluke Ohmmeter was used for measurement of the conductive contact time. The membrane size was cut to 0.75 in. square, and placed in the test holder, followed by addition of 400 grams of weight on the cover 104 for holding the samples in place. The opening 106 in cover 104 was 3/8" in diameter and used to make the tests shown in Table 2. Readings were made after 1 minute of contact with a drop of surfactant.

In this example, the stack of polycarbonate membranes was used to determine the ability of different surfactant solutions listed in Table 2 to penetrate the surface of the polymer artificial skin models. Breakthrough time was determined as the time required for the surfactant solution to be transported through the membrane(s) and contact and bridge two pieces of 20 gage copper wire (wire gap about 0.0625 inch).

Table 2 shows the results for fifteen surfactant systems. (5% solutions based on the active ingredients) and their penetration capabilities on polycarbonate membrane stacks. A direct comparison was also made between the surfactant solutions and deionized water as shown. The 3ML-15638 emulsion was used at 80% activity, since it was considered to already be very dilute.

Table 2. Initial Solution Screening Tests Through Polycarbonate Membranes (0.75" diameter).

Sample #	Surfactant	% Active	Surfactant Breakthrough Time (sec.)
EEG-1	Actrosol C-75™ (clear)	5%	.01.
EEG-2	Actrosol C-85™ (disp.)	5%	.05
EEG-3	Actrosol OY-75™ (disp.)	5%	.01
EEG-4	Chembetaine C™ (disp.)	5%	.01
EEG-5	Desonic DA-4™ (disp.)	5%	300
EEG-6	DeSulf SLES 301™ (clear)	5%	230
EEG-7	Dodecyl Sulfate Sodium	5%	14

	Salt (clear)		
EEG-8	Fluorad FC-129™ (clear)	5%	.06
EEG-9	Lamepon S™ (clear)	5%	.02
EEG-10	Lodyne S-103A™ (clear)	5%	.01
EEG-11	Lodyne S-222N™ (disp.)	5%	.03
EEG-12	MayTein C™ (clear)	5%	19
EEG-13	MayTein CT (clear)	5%	.01
EEG-14	3ML-15638™ (emulsion (o/w))	80%	34
EEG-15	T-MAZ 20™ Peg Ester (clear)	5%	.02
EEG-16	Deionized Water	100%	154

Example 3

- After the tests in Example 2, there was a concern that the 3/8" diameter opening was too small and may allow wicking of the test material around the edges of the test material onto the copper wires. Therefore, a different cover was prepared having about a one (1) inch opening for additional tests as shown in Table 3. Results were about the same. This showed that that wicking was not a problem. However, the larger diameter opening was used in the remainder of the tests. In addition, larger test material samples were prepared that were about 1.25 inch in diameter. Except for the larger opening and larger test material sample size, test parameters were the same as in Example 2.

Table 3. 5% Surfactant Solution Penetration Tests on 5 Layers (Total thickness = 30 µm) of Polycarbonate Membranes

Sample #	Surfactant	% Active	Surfactant Breakthrough Time (seconds)
EEG-1 AC	Actrosol C-75™ (clear)	5%	0.01
EEG-3 AY	Actrosol OY-75™ (disp.)	5%	0.01
EEG-9 LS	Lamepon S™ (clear)	5%	0.02
EEG-12 MC	MayTein C™ (clear)	5%	14.67
EEG-13 MT	MayTein CT™ (clear)	5%	0.10
EEG-15 TM	T-Maz 20™ Peg Ester (clear)	5%	27.00

Example 4. (Nylon membrane skin model)

The membrane used as another polymeric model for human skin was a nylon Spectra/Mesh™ membrane having 8 µm openings with a thickness of 75 µm (Cat # 1465128). Three layers of the nylon membrane were used for a total thickness of 225 µm. The membranes were layered so that the visible directional lines apparent on the surface were oriented in different directions. One layer was oriented ↑, the lines in the second layer were oriented to the right →, and the lines in the third layer were at a 45 degree angle to these. This staggering of the layer weave was done to prevent direct lineup of any pores in the material and to make the material mimic the skin to a higher degree.

Test parameters and conditions were the same as in the previous Examples. The membrane size was cut to 1.25 in. diameter placed in test holder followed by addition of 400 grams of weight for hold down. The deionized water (EEG-16) in Table 2 required 1:32 minutes before the nylon surface even wetted.

Tests were preformed for the formulations shown in Table 4 with addition of a dye so that wicking area could be evaluated. This example definitely showed that with a larger 1" opening that wicking was not a problem. The largest wicking area was 0.81 inch in diameter for Sample EEG-3 AY.

Table 4. Surfactant Breakthrough Time VS Water Using Three Layers (225 µm) of Nylon as the Membrane Material.

Sample #	Surfactant	% Active	Surfactant Break-through Time (seconds)	Resistivity at 1 Minute Post Break-through (Ohm-meter)
EEG-1 AC	Actrosol C-75™ (clear)	5%	0.01	N/A
EEG-3 AY	Actrosol OY-75™ (disp.)	5%	>300	No contact

EEG-9 LS	Lamepon S™ (clear)	5%	0.01	222 K ohms
EEG-12 MC	MayTein C™ (clear)	5%	0.01	1.25 M ohms
EEG-13 MT	MayTein CT™ (clear)	5%	0.01	1.16 M ohms
EEG-15 TM	T-MAZ 20™ Peg Ester (clear)	5%	0.01	4.67 M ohms
EEG-16	Deionized Water	100%	>300	No contact

Example 5

Tests were performed as in Example 4 except an additional layer of the same nylon membrane was added. This time the weave of the nylon layers was staggered as follows for each layer (1)↑, (2)→, (3)↓, (4)←. Results are shown in Table 5. Sample EEG-3 AY was not tested because no contact had been made with only three layers of nylon. Results are shown in Table 5. No contact was noted for samples EEG-1 AC and EEG-15 TM. It was found that EEG-1 AC and EEG-15 TM penetrated three layers but no solution was found beyond the three layers. For comparison, all the solutions in Table 3 went through five layers of a Nuclepore polycarbonate membrane.

Table 5. Surfactant Break Through Time Vs Water (Total thickness 300 µm - Four Layer Nylon Membrane as Skin Model)

Sample #	Surfactant	% Active	Surfactant Break-through Time (seconds)	Surfactant solution Resistivity (Ohm-meter)
EEG-1 AC	Actrosol C-75™ (clear)	5%	>300	No contact
EEG-9 LS	Lamepon S™ (clear)	5%	0.01	0.9734 M ohms
EEG-12 MC	MayTein C™ (clear)	5%	0.01	2.37 M ohms
EEG-13 MT	MayTein CT™ (clear)	5%	0.01	5.12 M ohms
EEG-15 TM	T-MAZ™ Peg Ester (clear)	100%	>300	No contact

Actrosol C-75 (EEG-1 AC) penetrated through three layers, but in the fourth layer there was no solution observed. T-MAZ Peg Ester (EEG-15 TM)

penetrated through three layers, but in the fourth layer there was no solution observed.

Example 6 - Human Subject Test Results

- 5 The tests were performed with Nicolet AgCl electrodes (019-417600), a impedance meter, a digital storage scope and digital voltmeter. A reference electrode and ground electrode were place on the subject on the right side of the forehead and the shoulder respectively. Standard abrasion and 10/20 gel were used. The reference to ground measurements were verified throughout
- 10 the test. Impedance was 1 kOhm and the half-cell potential was 18-22mV. All tests were done on the forehead with no skin preparation. All test substances were placed via a swab onto the forehead followed by the AgCl electrode held in place with medical tape. A new site and electrode was used with each test.
- 15 All solutions used were standardized to 5 % activity. Maytein C™; Lamepon S™; Dodecyl Sulfate, Na salt; and Actrosol C75™ gave results better or equal to reference samples 10/10 gel, Ultrasonic gel, and Elefix™ gel.

Table 6. Example 6 Results

20

Sample/ Substance (5% Solutions)	Initial Impedance (kOhm)	90 sec. Impedance (kOhm)	Initial ½ cell Potential (mv)	90 sec ½ cell potential (mv)	Notes	Surfactant Solution Resistivity (Ohm- meter)	Surface Tension (dynes/ cm)
Dry	>550	>550			Failed		
29-14 TM/ T-Maz 20	58	80	120	120	Climbing	77.2	36.0
29-12 MT/ Maytein CT	29	72	83	83	Climbing	2.2	26.7
29-11 MC/ Maytein C	16	18	70	70	Stable	1.1	27.5
29-8 LS/ Lamepon S	16	25	5	13	Slight climb	127.8	32.5
29-7 DS/ Dodecyl sulfate Na	13	12.9	120	103	Stable	1.6	34.2
29-4 CC/ Chembetaine C	27	36	110	91	Stable	1.7	34.3
29-3 AY/ Actrosol OY 75	25	29.8	100	83	Stable	1.9	32.7

29-1 AC/ Actrosol C75	16	28	46	35	Slight climb	1.1	34
10/20 gel	18	25	40	40	Slight climb		
Saliva	16	>550	2		Failed		
Ultrasonic gel	33	34	109	106	Stable		
Elefix Gel	39	28	33	36	Declining		

Example 7 Further Human Subject Tests

The tests were performed with a Nicolet AgCl electrode (019-417600), impedance meter, digital storage scope and digital voltmeter. A reference electrode and ground electrode were placed on right side of the forehead and the shoulder respectively of a human subject. Standard skin abrasion and a 10/20 gel were used. The reference to ground measurements were verified throughout the test. Impedance was about 900 ohms and the half-cell potential was about 6-12 mV.

All tests of substances on the human subject were done on the forehead with no skin preparation. All test substances were injected using a syringe into a cup electrode secured to the forehead with medical tape. Second and third injections were done at five-minute intervals. A new site and electrode was used with each test. All Substances were tested. Testing began at 10:47 AM and was completed at 2:03 PM. The order shown below is not the order in which the substances were tested. All solutions used were standardized to 5 % activity.

Table 7. Example 7 - Results for Substance Impedance (k Ω)

Sample No. ⇒	29-1 AC	29-3 AY	29-4 CC	29-7 DS	29-8 LS	29-11 MC	29-12 MT (2)	29-14 TM
	Resist. (k Ω)	Resist (k Ω)	Resist (k Ω)	Resist (k Ω)	Resist (k Ω)	Resist (k Ω)	Resist (k Ω)	Resist (k Ω)
Dry 1 st	34	Off Scale	140	132	***	175	28	124
Injection Time (min)								
0	16	42	66.5	26	14 (1)	14	16	54
1	14	44	1305	22	13 (1)	12	22	41
2	12	54	11	20	13 (1)	14	25	39

3	12	59	10	20	13 (1)	14	29	39
4	11	60	10	21	9	15	32	37
5	10	64	10	20	10	16	35	36
2 nd Injection Time (min)								
0	10	26	9	15	5	8	9	33
1	9	27	8	15	6	9	10	31
2	8	29	8	15	6	9	10	30
3	8	31	8	15	6	10	12	30
4	7	32	8	15	6	10	13	30
5	7	35	8	14	6	11	13	30
3 rd Injection Time (min)								
0	6	21	6	11	4	7	6	32
1	6	25	7	12	5	8	6	32
2	6	27	7	11	6	9	7	33
Post Test (3)	Off Scale	84	--- (4)	11	8.4	55	13	88
½ cell Potential (mv)	60-67	108- 143	116-133	-102 to -133	-29 to +3	-86 to -92	-84 to -99	145-192
Test 1 Impedance (kOhm)	16-28	25-30	27-26	13	16-25	16-18	29-72	58-74
Test 2 Impedance (kOhm)	10-33			9-23	3-9	8-24		

(1) The initial measurement was between the wire lead under test and reference and was corrected when discovered. The 13 kOhm readings correlated to 9.4 kOhm when the wiring was corrected.

(2) The lead physical attachment was a problem initially. This substance is sensitive to movement of the wire lead.

(3) Post Test measurement was performed after the completion of all the tests. Substance 29-8 LS retained a low impedance almost one hour after its initial test. Substance 29-8 LS was a Lamepon S formulation from first lot (327961).

(4) Measurement not taken.

Substances 1-AC, 29-4 CC, and 29-8 LS gave the best results.

Substance 29-8 LS, Lamepon S, gave the best results. Substances 29-3 AY and 29-14 TM had the worst results and accordingly are deleted from further

consideration by this test.

Example 8 - Human Subject Test Results

The tests were performed with a Nicolet AgCl electrode (019-417600), impedance meter, digital storage scope and digital voltmeter. A reference electrode and ground electrode were placed on right side of the forehead and the shoulder respectively. Standard abrasion and 10/20 gel were used. The reference to ground measurements were verified throughout the test. Impedance was 900 ohms and the half-cell potential was 6-12 mV. All tests were done on the forehead with no skin preparation. All test substances were injected using a syringe into a cup electrode secured to the forehead with medical tape. Second and third injections were done at five-minute intervals. A new site and electrode was used with each test. Only substances 29-1 AY, 29-7 DS, 29-8 LS, and 29-11 MT, identified above were used in this test.

Table 8. Example 8 Results - Substance Impedance (kOhm)

Time (min) ↓	29-1 AC	29-7 DS	29-8 LS	29-11 MC
Dry	Off Scale	137	Off Scale	140
Initial Inject.				
0	30.4	33	40	32.0
1	30.7	23.3	10.3	24.5
2	--- (2)	22.2	10.3	22.5
3	30.6	22.0	10.3	21.4
4	32.2	21.8	10.3	21.7
5	32.8	21.5	10.5	22.7
Second Inject				
0	16.0	12.5	5	13.5
1	--- (2)	12.0	5.7	17.7
2	18.2	11.5	5.9	19.1
3	--- (2)	11.2	6.1	21.4
4	--- (2)	11.6	6.1	23.0
5	20.2	12.5	6.0	25.1
Third Inject				
0	10.0	8.6	3.2	7.5
1	13.8	13.7	3.9	7.8
2	14.4	14.2	4.1	11.2
Post Test (1)	70.0	--- (2)	6.1	--- (2)

½ Cell Potential (mv)	44-57	96-135	2-13	71-82
Previous Test Impedance (kOhm)	19-28	12-13	16-25	16-18

1 - Post Test measurement was performed after the completion of the next set of electrode tests. The measurement was generally made about 17 - 20 minutes after the first test.

2 - Measurement not taken.

5

The substances all improved with additional injections. The rising impedance seen in the last test was reduced. Although not wishing to be constrained by theory, this effect may be explained by evaporation. Substance 29-8 LS, which is Lamepon S™, showed the best results in terms of resistance. No adverse skin reactions were found with any of the substances.

10

Example 9

Tests for surfactant solution electrical conductivity were performed on additional solutions including solutions having different activities and various excipients. These solutions were then used for additional human subject tests described further below.

15

Excipients such as polyethylene glycol 200 (PEG 200), NaCl, citric acid, and combinations of PEG 200 and NaCl were added to some of the basic surfactant solutions. In addition, the surfactant solutions had varied amounts of activity as shown in Table 9 below. The important results shown in this table are that increasing concentrations (% activity) of Actrosol C75, Lamepon S™, MayTein C™, and MayTein CT™ resulted in lower resistivity in the surfactant solution.

20

25

Table 9. Surfactant Solution Conductivity Data

Date Prepared L g Bo k Ref. #	Formulation	Surfactant % Active	Excipient	Surfactant Solution (vwr) Conductivity mS/cm	Solution Resistivity Ohm-Meter
42-01 AC	Actrosol C75	20	None	38.2	0.2618
42-02 AC	Actrosol C75	30	None	60.5	0.1653
42-03 CC	Chembetaine C™	20	None	21.4	0.4673
42-04 CC	Chembetaine C™	30	None	12.9	0.7776
42-05 DS	Dodecyl Sulfate, Na salt	5	0.1% PEG 200	5.9	1.6978
42-06 DS	Dodecyl Sulfate, Na salt	5	0.01% NaCl	6.2	1.6129
42-07 DS	Dodecyl Sulfate, Na salt	5	0.01% Citric acid	6.0	1.6611
42-08 DS	Dodecyl Sulfate, Na salt	5	0.1% PEG 200, 0.01 % NaCl	6.2	1.6234
42-09 DS	Dodecyl Sulfate, Na salt	20	None	41.0	0.2439
42-10 DS	Dodecyl Sulfate, Na salt	30	None	66.2	0.1511
42-11 LS	Lamepon S™*	5	0.1% PEG 200	78.0	0.1282
42-12 LS	Lamepon S™*	5	0.01% NaCl	77.3	0.1294
42-13 LS	Lamepon S™*	5	0.01% Citric acid	76.4	0.1309
42-14 LS	Lamepon S™*	5	0.1% PEG 200, 0.01 % NaCl	77.4	0.1292
42-15 LS	Lamepon S™*	20	None	223.0	0.0448
42-16 LS	Lamepon S™* (viscous)	30	None	142.5	0.0702
42-17 MC	MayTein C™	20	None	45.0	0.2222
42-18 MC	MayTein C™	30	None	66.0	0.1515
42-19 MT	MayTein CT™	20	None	15.7	0.6382
42-20 MT	MayTein CT™	30	None	18.0	0.5571

Note: All Lamepon S-based formulations with an * were prepared from first lot number (327961)

5 Example 10

Test performed with Nicolet AgCl electrode (019-417600), impedance meter, and digital voltmeter. A reference electrode and ground electrode were placed on left side of the forehead and the shoulder of a human subject respectively. Standard abrasion and 10/20 gel were used. The reference to ground measurements were verified throughout the test. Impedance and the

half-cell potential of the reference electrode were recorded and are listed in the results below. All tests substances were injected using a syringe into a cup electrode secured to the forehead with medical tape. Second and third injections were done at five-minute intervals. A new site and electrode was
5 used with each test. The substances were tested in pairs. All Substances were tested. The impedance meter was wired to remain on as it was in previous tests.

The substances from Table 9 were used to obtain the data listed in Table 10.

10

Table 10. Example 10 Resistivity Results (kOhm - unless identified otherwise)

[illegible]

Table 10. (Continued) - Example 10 Resistivity Results
(kOhm - unless identified otherwise)

	42-1 AC	42-2 AC	42-3 CC	42-4 CC	42-17 MC	42-18 MC	42-5 DS *	42-6 DS *	42-7 DS **	42-8 DS **	42-9 DS *	42-10 DS *	42-11 LS	42-12 LS	42-13 LS	42-14 LS	42-15 LS *	42-16 LS *	42-19 MT	42-20 MT
3 rd Injection	3	4.5	5.5	3	4	3.5	4.2	2.6	---	---	0.9	1.8	0.7	0.9	0.9	2.7	---	---	2.1	3.4
1	4	5		3.2	4.8	4	4.5	2.8	---	---	1.7	1.7	1.1	1.3	1.3	3.8	---	---	2.6	3.4
2	3.8	6	6.6	3	4.8	4	4	2.8	---	---	1.7	1.7	1.2	1.3	1.4	3.8	---	---	2.9	3.2
3	4																2	1		
Post Test	14.4 (70 min)	12.7 (70 min)											9.8 (25 min)	1.7 (25 min)			2.7 (35 min)	0.9 (35 min)		
½ cell Potential (mv)	4 to 19	-20 to -7	65 to 72	39 to 62	22 to 34	-4 to 16	55 to 106	71 to 120	66 to 109	77 to 101	51 to 74	85 to 106	12 to 17	18 to 21	9 to 11	-1 to -22	-34 to -37	-32 to -46	20 to 36	14 to 29

Notes for Table 10

* Slight Skin Irritation

** Skin Irritation

(1) Holding the electrode lowered the impedance. Several of the substances were soapy in nature and loosened the tape quickly.

(2) Repositioned fingers

(3) No third injection was given because of irritation.

(4) No third injection was needed. (Good electrical conductivity that was not improving)

(5) The initial dry site impedance was 6 kOhm. A new site and electrode were selected. Although the electrode site was about two inches away and was above the initial solution test sites, a residual effect may have been the cause.

(6) Holding the electrode lowered the impedance.

(7) Loose electrode

(8) All substances were tested for impedance. Readings in the table are all in kOhm except the ½ cell potential which is in (mv).

The Lamepon S™ appeared to be the best substance. Additives did not enhance its performance.

The PRG 200 may reduce the performance of the Lamepon S™.

- 5 Increasing concentration did not necessarily improve the performance of all the surfactants, although the 30% concentration was better in the Chembetaine C™ and Lamepon S™.

Method for Examples 11 through 16

- 10 Additional solutions illustrated in Examples 11 to 16 below were tested as follows. Tests were performed with a Nicolet AgCl electrode (Part No. 019-417600), impedance meter, and digital voltmeter. A reference electrode and ground electrode were placed on the left side of a subject's forehead and the subjects shoulder respectively. Standard abrasion and 10/20 gel were used.
- 15 The reference to ground measurements were verified throughout the test. Impedance and the half-cell potential of the reference electrode were recorded and are listed in the results below. All tests of the substances were done on the subject's forehead with no skin preparation (e.g. no mechanical abrasion or chemical preparation). All test substances were swabbed on the
- 20 subjects head with a previously dry Q-tip. A cup electrode was secured to the subject's forehead with essentially non-conducting medical tape. A new forehead site and electrode were used with each test. The impedance meter was wired to remain on as it was in previous tests.

25 **Example 11**

Gelling agents were added to Lamepon S™ and one sample of Dodecyl Sulfate as shown in Table 11.

30

Table 11. - Gelling Agents

Substance Sample No.	Surfactant	Gelling Agent (% by weight)
A-LS	Lamepon S ^a	NONE
B-LS	Lamepon S ^b	NONE
C-LS	Lamepon S ^a	3.2% Silver Micro-spheres
D-LS	Lamepon S ^a	3.6% Carbepol EZ-2
E-LS	Lamepon S ^a	1.0% Carbepol EZ-2
F-LS	Lamepon S ^a	1.0% Cellosize QP300
G-LS	Lamepon S ^a	1.0% Natrosol 250 LR (HEC)
H-LS	Lamepon S ^a	0.5% PVP K90
I-LS	Lamepon S ^b 5% active	NONE
J-LS	Lamepon S ^a 5% active	0.1% PEG 200
K-LS	Lamepon S ^b 20% active	NONE
L-DS	Dodecyl Sulfate, Na salt 5%	0.1% PEG 200 0.01% NaCl

a - Lamepon STM solution prepared from second lot - 30.4% active

b - Lamepon STM solution prepared from first lot - 52.7% active.

5

Example 12

Table 12 illustrates impedance measurements as outlined above for test subject A1. The electrodes were attached to the subject as described above. The reference impedance was 1.5 to 1.8 kOhms.

10

Table 12. Impedance Measurements (Subject A)

Sample No.	Impedance (kOhm)	½ Cell Potential (mv)	½ Cell Standard Deviation
A-LS	8.6	- 8.04	7.02
B-LS	11.9	58.16	2.27
C-LS	12.9	- 4.20	6.29
D-LS	12.5	5.88	1.37
E-LS	12.6	0.86	4.04

Sample No.	Impedance (kOhm)	½ Cell Potential (mv)	½ Cell Standard Deviation
F-LS	9.6	- 1.86	6.51
G-LS	8.6	6.24	4.46
H-LS	17.6	4.58	2.62

The values for impedance reported in Table 12 represent the lowest values obtained nominally after 6-7 minutes of testing. The results from Table 12 show that the new lot of Lamepon S™ gelled with 1% Cellosize QP300 and 1% Natrosol 250LR (HEC). Samples F and G respectively gave impedance results similar to the new lot ungelled Lamepon S™ of Sample A. Thus gelation with these gelling agents did not appear to affect impedance. An additional observation from this table is that Samples A and B show that as the concentration of Lamepon increases above about 30% to about 53% the impedance begins to increase.

Example 13

Table 13 illustrates impedance measurements as outlined above for test subject B1. Half-cell potentials were measured and half-cell standard deviations calculated. The electrodes were attached to the subject as described above. The test compares ungelled solutions of surfactant with gelled solutions. The reference impedance was 1.5 to 2.1 kOhms.

Table 13. Impedance Measurements (Subject B1)

Sample No.	% Activity	Impedance (kOhm)	½ Cell Potential (mv)	½ Cell Standard Deviation
L-DS	5	9.2	92.54	4.70
F-LS	30.4	7.3	- 1.21	2.26
J-LS	5	8.0	1.0	3.68
G-LS	30.4	12.2	2.40	7.15

Tests reported in Table 13 were side by side tests of Sample L-DS which was a 5% active solution gelled with 0.1% PEG 200 and 0.1% NaCl, and

Sample F-LS material which was 30.4% active Lamepon S™ gelled with 1% Cellosize QP300. A second side by side test was also made for Sample J-LS which was a 5% active solution of Lamepon S™ gelled with 0.1% PEG 200 and Sample G-LS which was a 30.4% active solution of Lamepon S™ gelled with 1% Natrosol 250 LR (HEC).

The results also show that the gelled Lamepon S™ of high activity outperformed the gelled Dodecyl Sulfate of low activity. This is important because Samples 42-10DS and 42-9DS of Table 9 indicated improvement in impedance as the activity of dodecyl sulfate was reduced. The results for Table 13 thus indicate that Lamepon S™ will outperform dodecyl sulfate at both high and low concentrations.

Example 14

Table 14 illustrates additional impedance measurements as outlined above for test subject B1. Half-cell potentials were measured and half-cell standard deviations calculated. The electrodes were attached to the subject as described above. The test compares gelled solutions for several surfactant solutions. The reference impedance was 1.7 to 2.0 kOhms. This example illustrates impedance tests for several combinations of surfactants and gelling agents.

Results indicate that adding PEG as a gellant resulted in a low initial impedance that increased over time for both dodecyl sulfate and Lamepon S™ (Samples L-DS and J-LS. Lamepon S™ gelled with 1% Cellosize QP300 or 1% Natrosol 250 LR(HEC) did not exhibit this phenomena but remained stable. When the materials gelled with PEG were mixed with Lamepon S™ gelled with the Cellosize or the Natrosol the instability in the form of a steady rise in impedance was reduced to a slight rise. This indicates that PEG is not a preferred gellant while Natrosol and Cellosize are.

30

Table 14. Additional Impedance Measurements (Subject B1)

Sample No.	Impedance (kOhm)	½ Cell Potential (mv)	½ Cell Standard Deviation
L-DS	11*	94.34	12.71
L-DS & F-LS	5.7**	9.97	4.75
F-LS	5.9	7.40	3.28
J-LS	26.7*	-9.67	14.91
J-LS & G-LS	10.0**	4.57	10.50
G-LS	9.1	3.83	3.55

Reference: 1.7 to 2.0 kΩ

* Initial reading was the lowest impedance with a steady rise over the course of the test.

- 5 ** One minute reading was the lowest impedance with a slight rise over the course of the test.

Example 15

- Table 15 illustrates further impedance measurements as outlined above
- 10 for test subject B1. Half-cell potentials were measured and half-cell standard deviations calculated. The electrodes were attached to the subject as described above. The example illustrates three ungelled solutions of Lamepon S™ surfactant and one Lamepon S™ gelled solution. A minimum impedance was found at about 20% activity for ungelled Lamepon S™.
- 15 Addition of 1.0% Natrosol 250 LR (HEC) to the highest (30.4% active) Lamepon S™ resulted in reduced impedance. It would thus be generally expected that addition of this type of gelling agent will besides other effects also enhance impedance. The reference impedance was about 2.4 kOhms.

20 **Table 15. Further Impedance Measurements (Subject B1)**

Sample No.	% Activity	Impedance (kOhm)	½ Cell Potential (mV)	½ Cell Standard Deviation
I-LS	5	5.9*	18.61	5.93
K-LS	20	2.7*	-22.65	2.24
A-LS	30.4	8.7	3.47	4.63
G-LS	30.4	7.3	10.96	3.69

Reference 2.4 kΩ

* A second injection was done on the thin solutions at 10 minutes.

Example 16

Table 16 illustrates yet additional impedance measurements as outlined above for test subject C1. Half-cell potentials were measured and half-cell standard deviations calculated. The electrodes were attached to the subject as described above. The example was essentially a repeat of Example 15, except that a different test subject was used. As before a minimum impedance was again measured with Lamepon S™ at an activity level of about 20%. This time, however, the addition of 1.0% Natrosol 250 LR (HEC) to the Lamepon S™ with highest activity (30.4% active) resulted in increased impedance. This is perhaps best explained as an anomaly, since for this sample the ½ cell Standard Deviation was much larger than the ½ cell potential. The reference impedance was about 2.4 kOhms.

Table 16. Impedance Measurements (Subject C1)

Sample No.	% Activity	Impedance (kOhm)	½ Cell Potential (mV)	½ Cell Standard Deviation
I-LS	5	10.2	13.51	9.48
K-LS	20	3.6	-24.86	5.53
A-LS	30.4	4.8	3.43	5.91
G-LS	30.4	7.0	6.51	9.65

Reference 2.4 kΩ

Drug Delivery

Another embodiment of the invention includes the use of the surfactants disclosed herein to enhance the delivery of drug compounds through the skin or membranes. The surfactant acts as carrier for drugs having small molecules, drugs with large molecules including proteins and other biologics. The embodiment covers the promotion of small or large molecule drug compounds through the skin. The transport of the drugs could be aided by an electrical potential. Typically, the drugs are used alone or are attached to

carrier molecules that can be charged with an electrical bias so as to enhance the transport of the molecules.

The transport enhancing material can be used on or in a patch for enhanced delivery/transport of drugs or other compounds. As an additive to
5 topical creams or gels with drug compounds mixed into the gel.

Non-Invasive Diagnostics

Another embodiment of the invention includes the use of the materials herein as enhancers for non-invasive diagnostic methods. As has been discussed in
10 earlier sections above, the surfactants of the invention have been demonstrated to increase electrical conductivity across the skin. These surfactants are also beneficial in promoting the uptake, diffusion, and capture of interstitial fluids through the skin. The interstitial fluids can carry proteins, biochemicals and other molecules from the body that are indicative of a
15 disease state, physiological condition, or response to a pharmaceutical or other therapy.

As an example, the surfactants of the invention can be embodied in a band-aid that is attached to the skin and draws out interstitial fluid that carries out the analyte of interest. Such analytes may include, but are not
20 limited to, insulin, antibodies developed by the immune system in response to infectious diseases such as strep infections, anthrax, hepatitis, and so on, or biomolecules that indicate response to a drug such as clotting factors or absorption of the drug in appropriate amounts. The surfactant/band-aid type diagnostic system could include immuno-chromatographic assays that change
25 color and give a visible measurement to the patient or health care provider about the physiological system.

Figure 4 illustrates a noninvasive diagnostic patch 400 that can be used to monitor bodily functions. The patch 400 consists of holding material 401 that
30 may be adhesive 401A at the bottom thereof where it contacts the skin. A first layer 402 containing the materials disclosed in this invention is located to be in proximity to the skin. A second reactive layer 403 is located above the

first layer and contains chemicals that are reactive with materials extracted from the body through the skin. A third layer 404 contains indicators visible to a user. Enveloping protective layer 405 is optional and may be clear or have a clear window so that the indicators are visible to the user.

5 A patch similar to that disclosed in Figure 4 may be used to deliver chemicals, drugs or medicine to a patient. In that case, layer 404 is not needed and only a drug containing layer 403 and a layer 402 containing the materials according to the present invention are needed. In some
10 embodiments the chemicals, drugs or medications may be mixed the materials of the present invention, in which case only layer 402 would be needed.

Electrosurgery

 In electrosurgery, large patch electrodes are used to establish a return
15 path for monopolar surgery tools. If the patch electrodes do not make sufficiently good contact with the patient's skin, the surgery may be less effective and, much worse, the patient may receive severe burns at the return electrode from the heat generated due to the higher contact resistance at the return electrode.

20 Currently, the gels applied to electrosurgery return electrodes contain highly ionic compounds. Addition of a surfactant such as those identified herein should have at least two benefits. First, the surfactant will drive the ions through the skin faster than if the ionic compounds are used in the absence of the surfactant. This will permit use of the electrodes more quickly
25 after their application to the patient. The second benefit is that it should be possible to increase the amount of ions present in the skin layer, thus reducing the contact resistance between the electrode and the skin. This would further reduce the risk of burns and improve the performance of the surgical procedure. Important aspects of the invention herein include the use
30 of the surfactant to improve ionic penetration of the skin and the specific properties of the surfactant mixtures that optimize its performance specifically for skin.

Referring now to Figure 5, this figure illustrates a schematic side view of a typical electrosurgery patch 500. A metallic conductive material 501 typically provides the support for the patch 500. Around the outside of the bottom 503 of the patch, an adhesive strip 502 typically surrounds the patch
5 perimeter. Conductive material 504 typified by the conductive materials of the invention are placed inside this area. The conductive material may be in a matrix material or may be applied as a gel that adheres to the electrode conductor 501. Typically a connector 505 is found on the upper surface of the electrosurgery electrode for connection to the source of power by a wire
10 510.

While the forms of the invention herein disclosed constitute presently preferred embodiments, many others are possible. It is not intended herein to mention all of the possible equivalent forms or ramifications of the
15 invention. It is to be understood that the terms used herein are merely descriptive, rather than limiting, and that various changes may be made without departing from the spirit of the scope of the invention.